



Catechol 1,2-dioxygenase Antibody

Product Code	CSB-PA306913ZA01AWR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P83715
Immunogen	Recombinant Acinetobacter calcoaceticus Catechol 1,2-dioxygenase protein
Raised In	Rabbit
Species Reactivity	Acinetobacter calcoaceticus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Catechol 1,2-dioxygenase (EC 1.13.11.1) (Fragment)
Immunogen Species	Acinetobacter calcoaceticus
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.