



MSP-1 Antibody

Product Code	CSB-PA361416ZA01EWT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04932
Immunogen	Recombinant Plasmodium falciparum MSP-1 protein
Raised In	Rabbit
Species Reactivity	Plasmodium falciparum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Merozoite surface protein 1 (Merozoite surface antigens) (PMMSA) (p190), MSP-1
Immunogen Species	Plasmodium falciparum (isolate K1 / Thailand)
Lead Time	Made-to-order (12-14 weeks)
Target Names	MSP-1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.