



dnaK-A Antibody

Product Code	CSB-PA656399ZA01DZX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q37106
Immunogen	Recombinant Cyanophora paradoxa dnaK-A protein
Raised In	Rabbit
Species Reactivity	Cyanophora paradoxa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chaperone protein dnaK (HSP70) (Heat shock 70 kDa protein) (Heat shock protein 70), dnaK-A; dnaK-B
Immunogen Species	Cyanophora paradoxa
Lead Time	Made-to-order (12-14 weeks)
Target Names	dnaK-A
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.