



IFNA6 Antibody

Product Code	CSB-PA011042ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05013
Immunogen	Recombinant Human Interferon alpha-6 (IFNA6) (21-189AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon alpha-6 (IFN-alpha-6) (Interferon alpha-54) (Interferon alpha-K) (LeIF K), IFNA6
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Gene Names	IFNA6
Product Modify	Unmodified