



fap2 Antibody

Product Code	CSB-PA892561ZA01SXV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UTM9
Immunogen	Recombinant Schizosaccharomyces pombe fap2 protein
Raised In	Rabbit
Species Reactivity	Schizosaccharomyces pombe
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	L-saccharopine oxidase (EC 1.5.3.18), fap2
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Lead Time	Made-to-order (12-14 weeks)
Target Names	fap2
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.