



DDX53 Antibody

Product Code	CSB-PA773024ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TM3
Immunogen	Recombinant Homo sapiens DDX53 protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable ATP-dependent RNA helicase DDX53 (EC 3.6.4.13) (Cancer-associated gene protein) (Cancer/testis antigen 26) (CT26) (DEAD box protein 53) (DEAD box protein CAGE), DDX53, CAGE
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	DDX53
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.