



botB Antibody

Product Code	CSB-PA533025ZA01DUJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B1INP5
Immunogen	Recombinant Clostridium botulinum botB protein
Raised In	Rabbit
Species Reactivity	Clostridium botulinum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Botulinum neurotoxin type B (BoNT/B) (EC 3.4.24.69) (Bontoxilysin-B) [Cleaved into: Botulinum neurotoxin B light chain; Botulinum neurotoxin B heavy chain], botB
Immunogen Species	Clostridium botulinum (strain Okra / Type B1)
Lead Time	Made-to-order (12-14 weeks)
Target Names	botB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.