



ruvC Antibody

Product Code	CSB-PA453002ZA01AZF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B2UP63
Immunogen	Recombinant Akkermansia muciniphila ruvC protein
Raised In	Rabbit
Species Reactivity	Akkermansia muciniphila
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Crossover junction endodeoxyribonuclease RuvC (EC 3.1.22.4) (Holliday junction nuclease RuvC) (Holliday junction resolvase RuvC), ruvC
Immunogen Species	Akkermansia muciniphila (strain ATCC BAA-835 / Muc)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	ruvC
Product Modify	Unmodified