



PPP1CA Antibody

Product Code	CSB-PA840958ZA01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WMS6
Immunogen	Recombinant Canis lupus familiaris PPP1CA protein
Raised In	Rabbit
Species Reactivity	Canis lupus familiaris
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase PP1-alpha catalytic subunit (PP-1A) (EC 3.1.3.16), PPP1CA
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Lead Time	Made-to-order (12-14 weeks)
Target Names	PPP1CA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.