




porA Antibody

Product Code	CSB-PA320687ZA01NGG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DH58
Immunogen	Recombinant Neisseria meningitidis serogroup B porA protein (20-392aa)
Raised In	Rabbit
Species Reactivity	Neisseria meningitidis serogroup B
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Neisseria meningitidis serogroup B (strain MC58)
Target Names	porA

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: porA Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 49 kDa

Observed band size: 49 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.