



porB Antibody

Product Code	CSB-PA333550ZA01NGG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30690
Immunogen	Recombinant Neisseria meningitidis serogroup B porB protein
Raised In	Rabbit
Species Reactivity	Neisseria meningitidis serogroup B
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major outer membrane protein P.IB (PIB) (Protein IB) (Class 3 protein) (Porin), porB
Immunogen Species	Neisseria meningitidis serogroup B (strain MC58)
Lead Time	Made-to-order (12-14 weeks)
Target Names	porB
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.