



nadA Antibody

Product Code	CSB-PA375490ZA01NEX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A1KVN9
Immunogen	Recombinant Neisseria meningitidis serogroup C / serotype 2a nadA protein
Raised In	Rabbit
Species Reactivity	Neisseria meningitidis serogroup C / serotype 2a
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Quinolinate synthase A (EC 2.5.1.72), nadA
Immunogen Species	Neisseria meningitidis serogroup C / serotype 2a (strain ATCC 700532 / DSM 15464 / FAM18)
Lead Time	Made-to-order (12-14 weeks)
Target Names	nadA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.