



walk Antibody

Product Code	CSB-PA409046ZA01FLG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6QD58
Immunogen	Recombinant Staphylococcus aureus walk protein
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sensor protein kinase Walk (EC 2.7.13.3), walk, vick
Immunogen Species	Staphylococcus aureus (strain Newman)
Lead Time	Made-to-order (12-14 weeks)
Target Names	walk
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.