



dltA Antibody

Product Code	CSB-PA336594ZA01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P39581
Immunogen	Recombinant Bacillus subtilis dltA protein
Raised In	Rabbit
Species Reactivity	Bacillus subtilis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	D-alanine--D-alanyl carrier protein ligase (DCL) (EC 6.2.1.-) (D-alanine--poly(phosphoribitol) ligase subunit 1) (D-alanine-activating enzyme) (DAE), dltA, dae
Immunogen Species	Bacillus subtilis (strain 168)
Lead Time	Made-to-order (12-14 weeks)
Target Names	dltA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.