



# enoA Antibody

<b>Product Code</b>	CSB-PA857046ZA01NGS
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96X30
<b>Immunogen</b>	Recombinant Neosartorya fumigata enoA protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Neosartorya fumigata
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Enolase (EC 4.2.1.11) (2-phospho-D-glycerate hydro-lyase) (2-phosphoglycerate dehydratase) (allergen Asp f 22), enoA
<b>Immunogen Species</b>	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Gene Names</b>	enoA
<b>Product Modify</b>	Unmodified