




H Antibody

Product Code	CSB-PA357366ZA01MCQ1
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08362
Immunogen	Recombinant Measles virus H protein (303-617aa)
Raised In	Rabbit
Species Reactivity	Measles virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Measles virus (strain Edmonston) (MeV) (Subacute sclerose panencephalitis virus)
Target Names	H

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes:H Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 39.4 kDa

Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.