




Abrin-a Antibody

Product Code	CSB-PA319970ZA01AAC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11140
Immunogen	Recombinant Abrus precatorius Abrin-a protein (1-251aa)
Species Reactivity	Abrus precatorius
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Abrus precatorius (Indian licorice) (Glycine abrus)
Research Area	Others
Target Names	Abrin-a

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: Abrin-a antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 33 kDa

Observed band size: 33 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.