



N Antibody

Product Code	CSB-PA494341ZA01MQM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B8PZP3
Immunogen	Recombinant Measles virus N protein
Raised In	Rabbit
Species Reactivity	Measles virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Nucleocapsid protein) (NP) (Protein N), N, NP
Immunogen Species	Measles virus (strain Edmonston-Schwarz vaccine) (MeV) (Subacute sclerose panencephalitis virus)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	N
Product Modify	Unmodified