



Oxalate oxidase 1 Antibody

Product Code	CSB-PA332739ZA01HWQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45850
Immunogen	Recombinant Hordeum vulgare Oxalate oxidase 1 protein
Raised In	Rabbit
Species Reactivity	Hordeum vulgare
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxalate oxidase 1 (EC 1.2.3.4) (Germin)
Immunogen Species	Hordeum vulgare (Barley)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified