



# BB3856 Antibody

<b>Product Code</b>	CSB-PA358313ZA01BTY
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0A321
<b>Immunogen</b>	Recombinant Bordetella bronchiseptica BB3856 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Bordetella bronchiseptica
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Azurin
<b>Immunogen Species</b>	Bordetella bronchiseptica (strain ATCC BAA-588 / NCTC 13252 / RB50) (Alcaligenes bronchisepticus)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	BB3856
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.