



P Antibody

Product Code	CSB-PA319601ZA01HGG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13846
Immunogen	Recombinant Heron hepatitis B virus P protein
Raised In	Rabbit
Species Reactivity	Heron hepatitis B virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein P [Includes: DNA-directed DNA polymerase (EC 2.7.7.7); RNA-directed DNA polymerase (EC 2.7.7.49); Ribonuclease H (EC 3.1.26.4)], P
Immunogen Species	Heron hepatitis B virus (HHBV)
Lead Time	Made-to-order (12-14 weeks)
Target Names	P
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.