



ret Antibody

Product Code	CSB-PA303691ZA01APQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P67876
Immunogen	Recombinant Aspergillus restrictus ret protein
Raised In	Rabbit
Species Reactivity	Aspergillus restrictus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribonuclease mitogillin (EC 3.1.27.-) (Restrictocin), ret
Immunogen Species	Aspergillus restrictus
Lead Time	Made-to-order (12-14 weeks)
Target Names	ret
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.