



B2M Antibody

Product Code	CSB-PA002486ZA01SH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6QAT4
Immunogen	Recombinant Ovis aries B2M protein
Raised In	Rabbit
Species Reactivity	Ovis aries
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-2-microglobulin, B2M
Immunogen Species	Ovis aries (Sheep)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	B2M
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.