



Minor capsid protein VP2 Antibody

Product Code	CSB-PA360956ZA01BGY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03094
Immunogen	Recombinant BK polyomavirus Minor capsid protein VP2 (2-351AA)
Raised In	Rabbit
Species Reactivity	BK polyomavirus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Minor capsid protein VP2 (Minor structural protein VP2)
Immunogen Species	BK polyomavirus (BKPyV) (Human polyomavirus 1)
Lead Time	Made-to-order
Product Modify	Unmodified