



dnapkcs Antibody

Product Code	CSB-PA703860ZA01DKK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q54UC0
Immunogen	Recombinant Dictyostelium discoideum dnapkcs protein
Raised In	Rabbit
Species Reactivity	Dictyostelium discoideum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA-dependent protein kinase catalytic subunit (DNA-PK catalytic subunit) (DNA-PKcs) (EC 2.7.11.1), dnapkcs
Immunogen Species	Dictyostelium discoideum (Slime mold)
Lead Time	Made-to-order (12-14 weeks)
Target Names	dnapkcs
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.