



UGT89C1 Antibody

Product Code	CSB-PA873303ZA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9LNE6
Immunogen	Recombinant Arabidopsis thaliana UGT89C1 protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-glycosyltransferase 89C1 (EC 2.4.1.-) (Flavonol 7-O-rhamnosyltransferase) (UDP-rhamnose: flavonol 7-O-rhamnosyltransferase), UGT89C1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order (12-14 weeks)
Target Names	UGT89C1
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.