



Dlat Antibody

Product Code	CSB-PA804374ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8BMF4
Immunogen	Recombinant Mus musculus Dlat protein
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dihydrolipoyllysine-residue acetyltransferase component of pyruvate dehydrogenase complex, mitochondrial (EC 2.3.1.12) (Dihydrolipoamide acetyltransferase component of pyruvate dehydrogenase complex) (Pyruvate dehydrogenase complex component E2) (PDC-E2) (PDCE2), Dlat
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (12-14 weeks)
Target Names	Dlat
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.