



# UGT71C1 Antibody

<b>Product Code</b>	CSB-PA526868ZA01DOA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O82381
<b>Immunogen</b>	Recombinant Arabidopsis thaliana UGT71C1 protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Arabidopsis thaliana
<b>Tested Applications</b>	ELISA, WB (ensure identification of antigen)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	UDP-glycosyltransferase 71C1 (EC 2.4.1.-) (Flavonol 3-O-glucosyltransferase UGT71C1) (EC 2.4.1.91) (Flavonol 7-O-glucosyltransferase UGT71C1), UGT71C1
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Lead Time</b>	Made-to-order (12-14 weeks)
<b>Target Names</b>	UGT71C1
<b>Product Modify</b>	Unmodified
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.