



UGT78D1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA871008ZA01DOA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9S9P6 |
| Immunogen | Recombinant Arabidopsis thaliana UGT78D1 protein |
| Raised In | Rabbit |
| Species Reactivity | Arabidopsis thaliana |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | UDP-glycosyltransferase 78D1 (EC 2.4.1.-) (Flavonol-3-O-glucoside L-rhamnosyltransferase) (EC 2.4.1.159) (UDP-rhamnose:flavonol 3-O-glucoside rhamnosyltransferase), UGT78D1 |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Lead Time | Made-to-order (12-14 weeks) |
| Target Names | UGT78D1 |
| Product Modify | Unmodified |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |