



traY Antibody

Product Code	CSB-PA327259ZA01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P33789
Immunogen	Recombinant Escherichia coli traY protein
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Relaxosome protein TraY (Fragment), traY
Immunogen Species	Escherichia coli
Lead Time	Made-to-order (12-14 weeks)
Target Names	traY
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.