



SREK1 Antibody

| | |
|---------------------|---|
| Product Code | CSB-PA837871ZA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8WXA9 |
| Immunogen | Recombinant Homo sapiens SREK1 protein |
| Raised In | Rabbit |
| Species Reactivity | Homo sapiens |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Splicing regulatory glutamine/lysine-rich protein 1 (Serine/arginine-rich-splicing regulatory protein 86) (SRrp86) (Splicing factor, arginine/serine-rich 12) (Splicing regulatory protein 508) (SRrp508), SREK1, SFRS12 SRRP86 |
| Immunogen Species | Homo sapiens (Human) |
| Lead Time | Made-to-order (12-14 weeks) |
| Gene Names | SREK1 |
| Product Modify | Unmodified |