




ftnA Antibody

Product Code	CSB-PA009053ZA01HUV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52093
Immunogen	Recombinant Helicobacter pylori Bacterial non-heme ferritin (ftnA) (1-167AA)
Species Reactivity	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
Target Names	ftnA

Image



Western Blot

Positive WB detected in: CSB-EP009053HUV(20ng), CSB-EP009053HUV(10ng)

All lanes: ftnA antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 24.8kDa

Observed band size: 24.8kDa

Exposure time: 5s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.