



Pollen allergen Lol p 1 Antibody

Product Code	CSB-PA321329ZA01LQK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14946
Immunogen	Recombinant Lolium perenne Pollen allergen Lol p 1 protein
Raised In	Rabbit
Species Reactivity	Lolium perenne
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pollen allergen Lol p 1 (Allergen Lol p I) (Allergen R7) (allergen Lol p 1)
Immunogen Species	Lolium perenne (Perennial ryegrass)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.