



RPS13 Antibody

Product Code	CSB-PA816848ZA01DTL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8I7D6
Immunogen	Recombinant Ciona intestinalis RPS13 protein
Raised In	Rabbit
Species Reactivity	Ciona intestinalis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	40S ribosomal protein S13, RPS13
Immunogen Species	Ciona intestinalis (Transparent sea squirt) (Ascidia intestinalis)
Lead Time	Made-to-order (12-14 weeks)
Target Names	RPS13
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.