




Htra1 Antibody

Product Code	CSB-PA879501ZA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9R118
Immunogen	Recombinant Mouse Serine protease HTRA1 (Htra1), partial (141-480AA)
Raised In	Rabbit
Species Reactivity	Mus musculus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease HTRA1 (EC 3.4.21.-) (High-temperature requirement A serine peptidase 1) (Serine protease 11), Htra1, Htra Prss11
Immunogen Species	Mus musculus (Mouse)
Target Names	Htra1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: Htra1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 39 kDa

Observed band size: 40 kDa

Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.