



G Antibody

Product Code	CSB-PA323293ZA01HPX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20896
Immunogen	Recombinant Human respiratory syncytial virus B G protein
Raised In	Rabbit
Species Reactivity	Human respiratory syncytial virus B
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major surface glycoprotein G (Attachment glycoprotein G) (Membrane-bound glycoprotein) (mG), G
Immunogen Species	Human respiratory syncytial virus B (strain 18537)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	G
Product Modify	Unmodified