




EBNA3 Antibody

Product Code	CSB-PA319032ZA01EFA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12977
Immunogen	Recombinant Epstein-Barr virus EBNA3 protein (1-138aa)
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Epstein-Barr virus
Target Names	EBNA3

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: EBNA3 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 22.2 kDa

Observed band size: 33 kDa