



SH3PXD2A Antibody

Product Code	CSB-PA735983ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TCZ1
Immunogen	Recombinant Homo sapiens SH3PXD2A protein
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 and PX domain-containing protein 2A (Adapter protein TKS5) (Five SH3 domain-containing protein) (SH3 multiple domains protein 1) (Tyrosine kinase substrate with five SH3 domains), SH3PXD2A, FISH KIAA0418 SH3MD1 TKS5
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order (12-14 weeks)
Target Names	SH3PXD2A
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.