



hykk Antibody

Product Code	CSB-PA428180ZA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8WFT6
Immunogen	Recombinant Danio rerio hykk protein
Raised In	Rabbit
Species Reactivity	Danio rerio
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydroxylysine kinase (5-hydroxy-L-lysine kinase) (EC 2.7.1.81), hykk, agphd1
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	hykk
Product Modify	Unmodified