



SAA1 Antibody

Product Code	CSB-PA020656ZA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35541
Immunogen	Recombinant Bos taurus SAA1 protein
Raised In	Rabbit
Species Reactivity	Bos taurus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum amyloid A protein (SAA) [Cleaved into: Amyloid protein A (Amyloid fibril protein AA)], SAA1, SAA
Immunogen Species	Bos taurus (Bovine)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	SAA1
Product Modify	Unmodified