



Genome polyprotein Antibody

Product Code	CSB-PA27089ZA01TJM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ICI2
Immunogen	Recombinant Turnip mosaic virus Genome polyprotein protein
Raised In	Rabbit
Species Reactivity	Turnip mosaic virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polyprotein
Immunogen Species	Turnip mosaic virus (strain Japanese) (TuMV)
Lead Time	Made-to-order (12-14 weeks)
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.