



BRLF1 Antibody

Product Code	CSB-PA662061ZA01EFC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3KSS7
Immunogen	Recombinant Epstein-Barr virus BRLF1 protein
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Replication and transcription activator (Rta) (Immediate-early protein Rta)
Immunogen Species	Epstein-Barr virus (strain GD1) (HHV-4) (Human herpesvirus 4)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	BRLF1
Product Modify	Unmodified