



S100a4 Antibody

Product Code	CSB-PA020632ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05942
Immunogen	Recombinant Rat Protein S100-A4 (S100a4) (2-101AA)
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A4 (Metastasin) (Nerve growth factor-induced protein 42A) (P9K) (Placental calcium-binding protein) (S100 calcium-binding protein A4), S100a4
Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order
Gene Names	S100a4
Product Modify	Unmodified