



S100a4 Antibody

Product Code	CSB-PA020632ZA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05942
Immunogen	Recombinant Rattus norvegicus S100a4 protein
Raised In	Rabbit
Species Reactivity	Rattus norvegicus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A4 (Metastasin) (Nerve growth factor-induced protein 42A) (P9K) (Placental calcium-binding protein) (S100 calcium-binding protein A4), S100a4
Immunogen Species	Rattus norvegicus (Rat)
Lead Time	Made-to-order (12-14 weeks)
Gene Names	S100a4
Product Modify	Unmodified