



yst Antibody

Product Code	CSB-PA339176ZA01YAE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31518
Immunogen	Recombinant Yersinia kristensenii yst protein
Raised In	Rabbit
Species Reactivity	Yersinia kristensenii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat-stable enterotoxin, yst
Immunogen Species	Yersinia kristensenii
Lead Time	Made-to-order (12-14 weeks)
Target Names	yst
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.