



acm Antibody

Product Code	CSB-PA326536ZA01SOH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25310
Immunogen	Recombinant Streptomyces globisporus acm protein
Raised In	Rabbit
Species Reactivity	Streptomyces globisporus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysozyme M1 (EC 3.2.1.17) (1,4-beta-N-acetylmuramidase M1), acm
Immunogen Species	Streptomyces globisporus
Lead Time	Made-to-order (12-14 weeks)
Target Names	acm
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.