



RHO Antibody

Product Code	CSB-PA019681ZA01PI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O18766
Immunogen	Recombinant Pig Rhodopsin (RHO), partial (1-36AA)
Raised In	Rabbit
Species Reactivity	Sus scrofa
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rhodopsin, RHO, RHO1
Immunogen Species	Sus scrofa (Pig)
Lead Time	Made-to-order
Gene Names	RHO
Product Modify	Unmodified