



VACWR088 Antibody

Product Code	CSB-PA362178ZA01VAI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07612
Immunogen	Recombinant Vaccinia virus VACWR088 protein
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Alias	Protein L1 (Virion membrane protein M25)
Immunogen Species	Vaccinia virus (strain Western Reserve) (VACV) (Vaccinia virus (strain WR))
Lead Time	Made-to-order (12-14 weeks)
Gene Names	VACWR088
Product Modify	Unmodified