



LOXL1 Antibody

Product Code	CSB-PA013040ZA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08397
Immunogen	Recombinant Human Lysyl oxidase homolog 1 (LOXL1) (95-574AA)
Raised In	Rabbit
Species Reactivity	Homo sapiens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysyl oxidase homolog 1 (EC 1.4.3.-) (Lysyl oxidase-like protein 1) (LOL), LOXL1, LOXL
Immunogen Species	Homo sapiens (Human)
Lead Time	Made-to-order
Gene Names	LOXL1
Product Modify	Unmodified