



disA Antibody

Product Code	CSB-PA762077ZA01MLA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q743W9
Immunogen	Recombinant Mycobacterium paratuberculosis disA protein
Raised In	Rabbit
Species Reactivity	Mycobacterium paratuberculosis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA integrity scanning protein DisA (Cyclic di-AMP synthase) (c-di-AMP synthase) (Diadenylate cyclase) (EC 2.7.7.85), disA
Immunogen Species	Mycobacterium paratuberculosis (strain ATCC BAA-968 / K-10)
Lead Time	Made-to-order (12-14 weeks)
Target Names	disA
Product Modify	Unmodified
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.